



Detector Developments, Strategies and Perspectives At ESRF

Heinz Graafsma

People who do the work:

- Detectors Pool: Menhard Kocsis, Marc Diot
- Special Detectors: John Morse, Thierry Martin, Cyril Ponchut, Menhard Kocsis
- Analogue Electronics: Jean-Claude Labiche
- Digital Electronics: Christian Herve
- Software group
- Mechanical engineering group
- Beamline scientists

Topics

- Structure at the ESRF
- Focus points
- Strategies
- High energy and high spatial resolution
- Avalanche Photo Diodes (APD's)
- Gas filled detectors
- New 2 D systems

Structure at the ESRF

- Detector Pool: off the shelf
- Special Detectors: consultation, tests, small developments
- Analogue/Transient elec: Frelon CCD
- Control elec: Specific electronics modules
- Digital elec (CS): Gas filled, drift diodes

Focus points

- High spatial resolution imaging (phosphors)
- Sub-milliseconds imaging (GFD's)
- New 2D systems (a-Si, pixel detectors, CMOS imagers, etc.)
- Beam monitoring
- High count rate 0D detectors (APD's, YAP:Ce).

Strategies

- Separate loan service from development.
- Support comes first, Developments second.
- Concentrate on what has most impact.
- Buy if you can, assemble if you have to, develop if you are allowed to.
- Collaborate: other SR, HEP, Space, Industry,...

HIGH ENERGIES IMAGING PERFORMANCE




Image depth:

14 bit

Pixel size:

14 x 14 μm²

Input pixel size:

0.47 μm or 1.4 μm

Field of view:

1x1 mm² or 2.9x2.9 mm²

Spatial resolution:

4 μm @ 65keV

Typical acquisition time with kinetic pipeline mode

	Binning	CCD+Exposure(45ms)	Time for 1000 frames
2048 x 128	no	13.6 fps	73 sec
2048 x 256	no	9.7 fps	104 sec
2048 x 256	2x2	13.3 fps	75 sec

DETECTOR on ID15



AVALANCHE PHOTODIODE

Design and construction of a fast counter for X-rays.

Made with Silicon Avalanche PhotoDiode.

- Energy range : $3 \text{ keV} < E_{\text{X-ray}} < 30 \text{ keV}$ (limited by thickness)
- Counting rate: $\sim 10^7 \text{ cps}$
- Time resolution: $\sim 1 \text{ ns}$
- Dark noise: $\sim 0.01 \text{ cps}$
- Energy resolution: $\sim 20 \% @ 24 \text{ keV}$
 $\sim 39 \% @ 12 \text{ keV}$
- Single control module, simple user interface

AVALANCHE PHOTODIODE



Head = APD + Pre-amplifier



7 Heads of detector available

- Hamamatsu
 - Two 5x3mm² 135 μ m available
 - ϕ =3mm 135 μ m (proto)
- EGG
 - Five 5x5mm² 110 μ m available
 - 10x10mm² 110 μ m (future)


Acquisition system : ACE (APD Controller Electronic)

- Principle of use: amplitude (mV) \Leftrightarrow energy(eV)
 - 1 counter, 2 thresholds (high and low) for level discrimination
 - Counter with low level only = integral counter.
 - Counter with low-high level = counter in energy range.

GAS-FILLED DETECTORS

- Activities in 2002
 - Development of 1D 10 x 200 mm² detector.
 - Micro-ionization chambers.
 - Further development of GEM technology (with CERN)
- Ongoing projects
 - Collaboration in development of Parallel DAQ based 2D detectors.
 - Position sensitive ionization chambers for beam position monitoring.
 - Obtained μm sensitivity with quadrupole chamber

Position sensitive quadrupole ionization chamber




PIXEL DETECTORS: MEDIPIX

Medipix is a collaborative project managed by CERN and involving several European institutes

The Medipix collaboration develops single photon counting area detector prototypes and evaluates them in :

- medical imaging
- non-destructive testing
- materials science
- nuclear decommissioning / gamma imaging

The ESRF joined the Medipix collaboration in August 2000



64x64 square pixels, 170 μm pitch
1.18 cm^2 sensitive area

PIXEL DETECTORS: MEDIPIX

SAXS: ID10

Beam : 8.12 keV, 20x20 μm


Sample : PMMA colloidal suspension

Acquisition : 100 x 0.2 sec exposures

Processing : dark subtraction

Maximum flux/pixel $\sim 10^5$ photons/s


6 orders of magnitude detected



PIXEL DETECTORS: MEDIPIX

SAXS : comparison Medipix/CCD/counter

Exposure times	
Medipix	20 s
Counter	30 min.
CCD	100 s



PIXEL DETECTORS

Ongoing Activities: IDEPHIX

Integrated project proposal for
European funding:

- Medical Imaging
- Non-destructive testing
- Safety inspection
- Synchrotron Radiation Science



CERN, ESRF, PSI


TEST OF CCD CAMERA DALSA: 1M60

Characteristics	Measured
Number of pixels	1024 x 1024
Pixel size	14 x 14 μm^2
Dark current	452 e ⁻ /pixel/sec @ 25°C
A/D converter	12 bit @ 4 x 20Mhz
Gain	50 e ⁻ /ADU
Dynamic Range	12 bit
Saturation	4095
Full Well capacity	204750 e ⁻
Readout noise	0.96 ADU
Readout time	60 fps
Storage	13sec → 800 frames




ANALOGUE AND TRANSIENT ELECTRONICS: CCD CAMERA

Scientific camera at ESRF




AVALANCHE PHOTODIODE

Real device
“Reach-Through” APD




Avalanche region Drift region




$APD \neq pin$  Internal Gain

AVALANCHE PHOTODIODE

Comparison of APD to Yttrium Aluminium Oxide Perovskite (YAP:Ce) scintillator cyberstar



@12.5 KeV Hamamatsu 5x3mm² 5.68ns
nonparalyzable model using log scales compare
to 1us NaI/PMT dead time.



	NaI(Tl)	YAP:Ce	APD
Availability	old	May 2003	May 2003
Active area	5cm ²	5cm ²	0.15cm ²
Energy range	5keV...	20...60keV	3...30keV
Energy resolution @ 22keV	28%	44%	30%
Dynamic with dead time correction	400Khz	2Mhz	50Mhz
Dead time	1μs	0.11μs	0.006 μs

ANALOGUE AND TRANSIENT ELECTRONICS: FRELON CAMERA

- Main achievements in 2002
 - Validation of concept for 2nd generation FRELON 16 bit (25000 grey levels), 40Mpixels/sec using DSP.
 - Design of a fast data acquisition board (2 Gbit/s serial data rate) dedicated to the 2nd generation FRELON. Possibility to extend to Dalsa camera.
 - In-house know-how: Gluing of one Taper on a CCD sensor (collaboration with J.Y. Massonnat Optics Group).
 - Kinetic Mode: 250fps for 16 lines: ID15, ID17.
- Ongoing Activities: FRELON 2nd generation
 - 2 prototypes for the end of 2003.
 - Ongoing purchase of 2 new tapers ~ 100 x 100 mm² (collaboration with ADSC).
- Camera on Beamline
 - 1k x 1k: ID11, ID19, probably ID22
 - 2k x 2k: BM5, ID2, ID13, ID17, ID19, ID22
 - Spares: two 1kx1k cameras, no 2kx2k camera

ANALOGUE AND TRANSIENT ELECTRONICS:

Comparison of FRELON to DALSA CAMERA

Characteristics	FRELON	DALSA
Number of pixels	2048 x 2048	1024 x 1024
Pixel size	14 x 14 μm^2	14 x 14 μm^2
Dark current	3 e ⁻ /pixel/sec @ -20°C	452 e ⁻ /pixel/sec @ 25°C
A/D converter	14 bit @ 4 x 5Mhz	12 bit @ 4 x 20Mhz
Gain	20 e ⁻ /ADU	50 e ⁻ /ADU
Dynamic Range	14 bit	12 bit
Saturation	16383	4095
Full Well capacity	320000 e ⁻	204750 e ⁻
Readout noise	1 ADU	0.96 ADU
Readout time	5 fps	60 fps
Storage		13sec → 800 frames

PIXEL DETECTORS: MEDIPIX

Ongoing Activities: Medipix-2

Medipix-1

Readout chip:

- 64x64 square pixels, 170 μm pitch
- Variable threshold
- 2 MHz/pixel count rate
- 15 bit counters
- Pixel threshold tuning

X-ray sensor:

- 1.18 cm^2 sensitive area

Medipix-2

Readout chip:

- 256 x 256 pixels, 55x55 μm pitch
- Energy windowing (2 thresholds)
- Positive or negative pulses (compatible with Si, CdTe, AsGa,...sensors)
- 1 MHz count rate
- 13 bit counters

X-ray sensor:

- 1.98 cm^2 sensitive area

Status

- First X-ray tests made at CERN in December 2002
- Tests planned at ESRF in 2003

PIXEL DETECTORS: MEDIPIX

^{55}Fe source (5.9 keV, 6.4 keV)


Vth steps : 5 mV

Exposure time : 10 s

Summing over all pixels

Low energy calibration mask

Energy resolution: possible



HIGH ENERGIES IMAGING


State-of-the-Art

FOLDED DETECTOR for Low Energy

- TOTAL MAGNIFICATION: $x2.3 \rightarrow x100$
- SCINTILLATOR: YAG:Ce^{3+} , LAG:Eu^{3+} or Tb^{3+}
 $1 \rightarrow 25 \mu\text{m}$
- SPATIAL RESOLUTION: $0.8 \mu\text{m}$

Luminescent screen
 $\text{Y}_3\text{Al}_5\text{O}_{12}:\text{Ce}^{3+}$ on
 $170 \mu\text{m}$ undoped YAG substrate

X-rays



Cooled CCD

Eyepiece


Tube lens

Mirror

Microscope objective




HIGH ENERGIES IMAGING



HIGH ENERGIES IMAGING

Ongoing activities

- New x10 reflecting objective (compromise between magnification and NA)
- Development of high-energy high-resolution imaging systems with 3 motorized zooms for ID11
- Thicker $\text{Lu}_2\text{O}_3:\text{Eu}^{3+}$ and $\text{Gd}_2\text{O}_3:\text{Eu}^{3+}$ scintillator for high energies (PLD)



TEST OF CCD CAMERA DALSA: 1M60

- **Advantages:**
 - High data rate 60fps
 - Low readout noise
 - Compact 94 x 94x 102 mm³
 - Good understanding of technical requirements by Dalsa
 - Camera for fast tomography
- **Disadvantages:**
 - Important dark current → short exposure time
 - Cooling system impossible
 - Quality of chip (5 columns defect)
 - 12 bit

