



NORTHERN ILLINOIS
UNIVERSITY



Materials Science Division
Argonne National Laboratory

The study of ferroelectric switching using x-ray synchrotron radiation

Carol Thompson

*Science with Microbeams
APS Scientific Advisory Cross-cut Review
January 21, 2004*

*This work was supported by the U. S. Department of Energy,
Basic Energy Sciences, under contract W-31-109-ENG-38, and
the State of Illinois, USA, under Higher Education Collaborative
Act.*

Outline

◆ Introduction

- ◆ *What is a ferroelectric*
- ◆ *Concentrate on epitaxial films*
- ◆ *Oxide perovskite system*

◆ Structural response of epitaxial ferroelectric to electric field

- ◆ *examples of dynamic studies*

◆ Summary and Conclusions

- ◆ *Domain studies, device studies, future studies require microbeams.*

◆ Collaborators



**NORTHERN ILLINOIS
UNIVERSITY**

- ◆ *Chris Gunderson (Physics, NIU)*
- ◆ *Marian Aanerud (Masters 2002, Physics, NIU)*






Materials Science Division
Argonne National Laboratory




- ◆ *Stephen Streitter (MSD, ANL)*
- ◆ *Brian Stephenson (MSD, ANL)*
- ◆ *G. -R Bai (MSD, ANL)*
- ◆ *W. K. Kee (XFD-XRP, ANL)*
- ◆ *Armon McPherson, (currently at Sandia)*

What are ferroelectrics?

- ◇ Spontaneous permanent electric polarization.
- ◇ Unit cell of crystal is non-centrosymmetric (charges separated)
- ◇ A macroscopic sample with net zero polarization
 - ◆ combination of microscopic polarized domains.



What are ferroelectrics?



Synchrotron techniques are well matched to the study of the ferroelectric systems

◆ Structure-property relationships control:

- ◆ dielectric, ferroelectric, piezoelectric, electrostrictive, pyroelectric and electro-optical properties*


for actuators, sensors, electro-optical switches, non-volatile memory elements, hi-K dielectric, detectors...

◆ Scattering and diffraction examine the structural aspects that control the properties


- ◆ Symmetry changes, orientation, lattice parameters, domains configurations*

Scattering example: fingerprints domain evolution

- ◆ Time-resolved scattering
 - ◆ 40 nm $\text{Pb}(\text{Ti,Zr})\text{O}_3$ film
 - ◆ 200 Hz





Difficult Easy





- ◆ Scattering profile can fingerprint the domain configuration in epitaxial films

High speed time-resolved Methods (BESSRC 12-ID-D)

- ◆ At each voltage, collect all scattering (area detector)
- ◆ Utilizes rocking curve of sample to “scan” q
 - ◆ Chopper synchronized (Hybrid fill: Singlet produces < 100 psec x-ray probe pulses)
 - ◆ Electrical stimulation of device synchronized/delayed so that sample is in particular electrical state during exposure






Close-up photograph of sample manipulation and contact region



*Spot size used: $5\mu\text{m} \times 5\mu\text{m}$
K-B mirror focus*


- ◆ *X-ray spot must be smaller than the device.*
- ◆ *And x-ray spot must be aligned with the device under electrical stimulation.*

More Pictures




Reciprocal Space Map 001



- ◆ Initial experiments: Focus on position of film Bragg peak region and its immediate neighborhood.



- ◆ Scattering shown for epitaxial films (thickness $\sim 250\text{nm}$) of PMN and PMN-PT



PMN₇-PT₃ Structural Response to a Step Voltage



- ◇ *Response (speed) limited by size of device, not by how fast we can measure with x-rays yet*
- ◇ *Smaller devices – smaller beams*

Summary and Conclusion

- ◇ *Structural techniques available at synchrotrons well suited to ferroelectric systems*
 - ◆ *And it's a growing field: see also other groups doing exciting studies of ferroelectric films and crystals using microdiffraction, x-ray topography, and reciprocal space mapping.*
- ◇ *Examples from our work:*
 - ◆ *Progress in development of techniques to study structural response at 100 psec time scale*
 - ◆ *Need to go to smaller devices, embedded devices*
 - Progress in switching studies: to 50 μm 'play' device: switching speed limited to $\sim 10\text{nsec}$*
 - Smaller devices allow faster switching*
 - Need for microbeam capabilities*


Preferred Domain Pinning

◆ Piezo-response atomic force microscopy:

- ◆ Recent direct observation of *preferred domain pinning* in fatigued ferroelectric films is reported using piezo-response atomic force microscopy.

Direct observation of inversely polarized frozen nanodomains in fatigued, ferroelectric memory capacitors, E. L. Colla, I. Stolichnov, P. E. Bradely, and N. Setter, **Appl Phys. Lett.** **82**, 1604 (2003).



Samples: Pt-PZT-Pt films.



phase amplitude
AFM- Piezoreponse image

Time-Resolved Synchrotron X-Ray Scattering


◆ Data taken on 250 nm thick PMN-PT film (PT ~ 30-35%)




Lattice response: time-resolved x-ray diffraction

◆ Lattice response on different time scales

- ◆ pulse with ~ 15 nsec rise time



- ◆ 6.3 kHz triangle



PMN₇-PT₃

